


# Technical Bulletin



## Improved Attachment and Appearance of Stop Dot on CS 53/63 Flex Strip

**S**onion is pleased to announce that we have improved the product quality of the CS 53/63 Flex Strip by improving the stop dot in relation to:

- Significant better adhesion of the stop dot to the flex strip
- More uniform shape and appearance of the stop dot.

The above improvements have been achieved by changing the dot material and improving the manufacturing process.

Qualification report is available and can be provided upon request.

The change will not have any influence on part numbers or ordering situation.

Should you, however, have any questions or need further information, please do not hesitate to contact your sales representative.

Dan Lundahl  
Quality Manager  
Medical Division