Data Sheet

Receiver 2023



Description

Miniature magnetic receiver (balanced armature type) for use in hearing aids with push-pull output stage.

Features


- Great BTE applications
- High output, maximum peak output 140 dB
- Balanced armature design

Mechanical data

Weight 0.94 gr. Case material Ni80Fe20

Solder pad material Sn96.5Ag3.0Cu0.5
Dimensions Refer to outline drawing

Product drawing - Dimensions in mm [inch]



Sonion reserves the right to make changes at any time to improve reliability, function or design, in order to provide the best product possible. Receivers series of this type can produce very high sound pressure levels. When such receivers are applied in hearing instruments or other communications equipment special attention should be paid to this capacity in order to prevent possible hearing damage.

DK: +45 4630 6666 USA: +1 952 543 8300 PRC: +86 512 6832 3401 NL: +31 20 6068 100



Data Sheet

Receiver 2023



Specifications
The acoustic termination consists of: 8 mm x 1 mm ID + 28 mm x 1.5 mm ID + 25 mm x 2 mm ID + 18 mm x 3 mm ID + 2 cc coupler. The electrical input is a 1.51 mA RMS constant current signal from a high impedance source, applied to the total winding. Environmental conditions: 23 °C (73.4F), 50 % RH.

B. 6 :				
Min	Тур	Max	Unit	Comments
110	112	114	dB	
111	113	115	dB	
113	115	117	dB	
750	875	1000	Hz	
125	127	129	dB	
1300	1600	1900	Hz	
116	118		dB	
1850	2000	2150	Hz	
125	127	129	dB	
2300	2550	2800	Hz	
116	119		dB	
2850	3000	3150	Hz	
121	123	125	dB	
3500	3700	3900	Hz	
109	112		dB	
4000	4150	4300	Hz	
113	116	119	dB	
4800	5100	5400	Hz	
97	100		dB	
5400	5700	5900	Hz	
105	110	115	dB	
		9	%	
		9	%	
130	140	143	dB	
	110 111 113 750 125 1300 116 1850 125 2300 116 2850 121 3500 109 4000 113 4800 97 5400 105	110 112 111 113 113 115 750 875 125 127 1300 1600 116 118 1850 2000 125 127 2300 2550 116 119 2850 3000 121 123 3500 3700 109 112 4000 4150 113 116 4800 5100 97 100 5400 5700 105 110	110 112 114 111 113 115 113 115 117 750 875 1000 125 127 129 1300 1600 1900 116 118 1850 2000 2150 125 127 129 2300 2550 2800 116 119 2850 3000 3150 121 123 125 3500 3700 3900 109 112 4000 4150 4300 113 116 119 4800 5100 5400 97 100 5400 5700 5900 105 110 115 9 9 9	110 112 114 dB 111 113 115 dB 113 115 117 dB 750 875 1000 Hz 125 127 129 dB 1300 1600 1900 Hz 116 118 dB 1850 2000 2150 Hz 125 127 129 dB 2300 2550 2800 Hz 116 119 dB 2850 3000 3150 Hz 121 123 125 dB 3500 3700 3900 Hz 109 112 dB 4000 4150 4300 Hz 113 116 119 dB 4800 5100 5400 Hz 97 100 dB 5400 5700 5900 Hz 105 110 1

Electric parameters	Min	Тур	Max	Unit	Comments
Impedance @ nominal	352	440	528	Ohm	Geometric avarage 1, 1.6 and 2.5 kHz
DC resistance @ 20°C	80	95	114	Ohm	

Additional parameters	Min	Тур	Max	Unit	Comments
Shock resistance	5500			g	90% survival rate with THD @ 1/2 peak frequency < 10%
Storage temperature range	-40		63	°C	

A positive voltage applied to the negative terminal (-) will result in an increase in pressure at the sound outlet.

Sonion reserves the right to make changes at any time to improve reliability, function or design, in order to provide the best product possible. Receivers series of this type can produce very high sound pressure levels. When such receivers are applied in hearing instruments or other communications equipment special attention should be paid to this capacity in order to prevent possible hearing damage.

DK: +45 4630 6666 USA: +1 952 543 8300 PRC: +86 512 6832 3401 NL: +31 20 6068 100




Data Sheet

Receiver 2023







Sonion reserves the right to make changes at any time to improve reliability, function or design, in order to provide the best product possible. Receivers series of this type can produce very high sound pressure levels. When such receivers are applied in hearing instruments or other communications equipment special attention should be paid to this capacity in order to prevent possible hearing damage.

DK: +45 4630 6666 USA: +1 952 543 8300 PRC: +86 512 6832 3401 NL: +31 20 6068 100

