


Technical Bulletin



PJ 84 Improved Marking Possibilities

To improve the color coding possibilities of the PJ 84 knob, Sonion has introduced a new Knob No. 2.

The mechanical dimensions are identical to Knob No. 1. However, as seen on the above illustration, this knob has three small cavities for color dot marking, located in each of the three areas which divide the knob.

These indentations allow for a blue/red color coding of a higher and more consistent quality compared to the previous Knob No. 1.

For this reason, we have also decided to discontinue Knob No. 1.

For current users of the PJ 84 with Knob No. 1, you will receive a new Sonion part number next time you order the PJ 84.

Morten Wagner
Product Line Manager
Medical Division