



Active Telecoils Color Coding Scheme

Introduction

The purpose of the color coding process is to improve product identification for all amplified telecoils and we have therefore added two color stripes on the hybrid encapsulate.

Description of Coding

Since Sonion offers a variety of amplified telecoils with potentially the same physical appearance, we have implemented a procedure to add two color stripes to the top of the black encapsulate. Currently, we have identified the use of six different colors that will best differentiate the telecoils.

Each amplified telecoil part number has been assigned a different color code. To protect the color marking from physical wear, we implement the color marking procedure before the top layer protective coating is applied.

We are providing this service to only our amplified telecoils and not to our passive coils. Passive coils have a variety of physical and electrical differences to help you more easily identify the part. As new amplified telecoil models are developed, we will continue to create new color code schemes to help you better identify the product.

Results

For the complete understanding to the color coding matrix, please see the below overview.

Amplified Telecoil Color Code Matrix		Color Marking	
Item no.	Type Name	Dot 1	Dot 2
9010-2000180	TA 20 AD 01	● Yellow	● Blue
9009-3005110	TA 20 AD 11	● Yellow	No line
9008-3009565	TA 20 AD 21	● Blue	● Blue
9010-2000341	TA 20 AE 01	● Yellow	● Yellow
9010-2002827	TA 20 AE 03	● Yellow	● Green
9008-3008312	TA 20 AE 21	● Blue	● Grey
9010-3002733	TA 20 CE 01	● Green	● Grey
9008-3016674	TA 32 CE 21	○ White	○ White
9010-3023812	TA 33 CE 01	○ White	● Yellow

Item no.	Type Name	Dot 1	Dot 2
9010-2000180	TA 20 AD 01	● Yellow	● Blue

Item no.	Type Name	Dot 1	Dot 2
9009-3005110	TA 20 AD 11	● Yellow	No line

Item no.	Type Name	Dot 1	Dot 2
9008-3009565	TA 20 AD 21	● Blue	● Blue

Item no.	Type Name	Dot 1	Dot 2
9010-2000341	TA 20 AE 01	● Yellow	● Yellow

Item no.	Type Name	Dot 1	Dot 2
9010-2002827	TA 20 AE 03	● Yellow	● Green

Item no.	Type Name	Dot 1	Dot 2
9008-3008312	TA 20 AE 21	● Blue	● Grey

Item no.	Type Name	Dot 1	Dot 2
9010-3002733	TA 20 CE 01	● Green	● Grey

Item no.	Type Name	Dot 1	Dot 2
9008-3016674	TA 32 CE 21	○ White	○ White

Item no.	Type Name	Dot 1	Dot 2
9010-3023812	TA 33 CE 01	○ White	● Yellow

Remark

Please contact your local sales representative for any questions you may have. We thank you for your business and are pleased to continually improve the quality and performance of our products.